



UNIVERSITY OF
CAMBRIDGE

Protein Purification Part 2

Day 4: Thursday 23rd March

Continuing from Yesterday

- In Part 1 we learnt about:
 - Lysis techniques
 - Affinity chromatography including affinity tags and resins
- In Part 2 we will learn about:
 - Other chromatography techniques
 - Size-exclusion
 - Ion exchange
 - HIC and others
 - Chromatography equipment
 - The dos and don'ts of good chromatography



Why might you need to do a 2-step purification?

- Isn't my protein pure enough after affinity chromatography?



Why might you need to do a 2-step purification?

- Isn't my protein pure enough after affinity chromatography?
- Possibly not, other protein contaminants can interfere with subsequent experiments
 - Tomorrow we'll discuss how to evaluate this




Why might you need to do a 2-step purification?

- Isn't my protein pure enough after affinity chromatography?
- Possibly not, other protein contaminants can interfere with subsequent experiments
 - Tomorrow we'll discuss how to evaluate this
- Also, there might be components in our buffer we need to remove
 - Although dialysis is one approach, 2-step purification is another




Size Exclusion Chromatography

- Separates based on size




Size Exclusion Chromatography

- Separates based on size
- Large proteins elute first
- Smaller ones are trapped in the resin and elute later




Size Exclusion Chromatography

- Separates based on size
- Large proteins elute first
- Smaller ones are trapped in the resin and elute later




Monitoring Size Exclusion Chromatography

- The chromatogram:




Monitoring Size Exclusion Chromatography

- The chromatogram:
 - Sample injection




Monitoring Size Exclusion Chromatography

- The chromatogram:
 - Sample injection
 - Void volume, V_0




Monitoring Size Exclusion Chromatography

- The chromatogram:
 - Sample injection
 - Void volume, V_0
 - Separation range
 - Large (aggregates) first




Monitoring Size Exclusion Chromatography

- The chromatogram:
 - Sample injection
 - Void volume, V_0
 - Separation range
 - Large (aggregates) first
 - Well-resolved proteins



Monitoring Size Exclusion Chromatography


- The chromatogram:
 - Sample injection
 - Void volume, V_0
 - Separation range
 - Large (aggregates) first
 - Well-resolved proteins
 - Small molecules



Monitoring Size Exclusion Chromatography

- The chromatogram:

- Note this is a simple “isocratic” flow
- One buffer, running continuously, no gradient, no elution buffer



Separation range of your column

- Different pore sizes to separate different proteins





Small molecules
able to
enter the pores



Separation range of your column

- Different pore sizes to separate different proteins


1. Ferritin	440 kDa
2. IgG	158 kDa
3. Albumin	66 kDa
4. Ovalbumin	44 kDa
5. Myoglobin	17 kDa




Separation range of your column

- Different pore sizes to separate different proteins

1. Ferritin	440 kDa
2. IgG	158 kDa
3. Albumin	66 kDa
4. Ovalbumin	44 kDa
5. Myoglobin	17 kDa



← Better separation of small proteins




← Better separation of large proteins

Separation range of your column

- For even bigger proteins
 - Superose 6 column

1. Thyroglobulin	660 kDa
2. Ferritin	440 kDa
3. Aldolase	160 kDa
4. Ovalbumin	44 kDa
5. Ribonuclease	14 kDa
6. Aprotinin	6.5 kDa



Types of Size Exclusion Chromatography

- Different column sizes for:
 - different amounts of protein
 - different resolution
 - different chromatography systems



Types of Size Exclusion Chromatography

- Different column sizes for:
 - different amounts of protein
 - different resolution
 - different chromatography systems



5/150



10/300



16/600



26/600



What if size doesn't separate my contaminants?

- Although SEC is excellent for removing large aggregates and exchanging small buffer components...
- ...sometimes the contaminants in our samples are the same size
- So we use different chromatography techniques to separate these





Ion Exchange Chromatography

- Separates based on charge





Ion Exchange Chromatography

- Separates based on charge





Ion Exchange Chromatography

- Separates based on charge
- Anion and cation exchange





Ion Exchange Chromatography

- Separates based on charge
- Anion and cation exchange
- How do I know what charge my protein is?



Ion Exchange Chromatography

- Separates based on charge
- Anion and cation exchange
- How do I know what charge my protein is?
- Using ProtParam as described on Day 2



My protein charge

- Protparam calculates a theoretical pI
- This is the pH where your protein is uncharged

Number of amino acids: 510

Molecular weight: 57364.86

Theoretical pI: 6.37



Amino acid composition:

Ala (A)	31	6.1%
Arg (R)	34	6.7%
Asn (N)	30	5.9%
Asp (D)	35	6.9%
Cys (C)	9	1.8%
Gln (Q)	24	4.7%
Glu (E)	27	5.3%
Gly (G)	32	6.3%
His (H)	9	1.8%
Ile (I)	23	4.5%
Leu (L)	45	8.8%
Lys (K)	25	4.9%
Met (M)	11	2.2%
Phe (F)	16	3.1%
Pro (P)	34	6.7%
Ser (S)	32	6.3%
Thr (T)	33	6.5%
Trp (W)	7	1.4%
Tyr (Y)	18	3.5%
Val (V)	35	6.9%
Pyl (O)	0	0.0%
Sec (U)	0	0.0%




My protein charge

- Protparam calculates a theoretical pI
- This is the pH where your protein is uncharged
- Above this it is -ve, below it is +ve

Number of amino acids: 510


Molecular weight: 57364.86

Theoretical pI: 6.37




My protein charge

- Protparam calculates a theoretical pI
- This is the pH where your protein is uncharged
- Above this it is -ve, below it is +ve




Ion exchange chromatography

- For this protein we could do:
 - Cation exchange in buffer @ pH 5
- OR
- Anion exchange in buffer @pH 8





Cation Exchange Chromatography

- Using cation exchange:
 - positively charged protein binds the column
 - negatively charged protein flows through





Cation Exchange Chromatography

- Using cation exchange:
 - positively charged protein binds the column
 - negatively charged protein flows through
 - Load the column in low [salt]





Cation Exchange Chromatography

- Using cation exchange:
 - positively charged protein binds the column
 - negatively charged protein flows through
 - Load the column in low [salt]
 - Negative proteins don't bind





Cation Exchange Chromatography

- Using cation exchange:
 - positively charged protein binds the column
 - negatively charged protein flows through
 - Load the column in low [salt]
 - Negative proteins don't bind
 - Elute using increasing [salt]



Cation Exchange Chromatography

- Using cation exchange:
 - positively charged protein binds the column
 - negatively charged protein flows through
 - Load the column in low [salt]
 - Negative proteins don't bind
 - Elute using increasing [salt]
 - Wash with high [salt]



Other chromatography techniques


- Hydrophobic interaction chromatography



Other chromatography techniques


- Hydrophobic interaction chromatography

- Separates based on hydrophobicity



Other chromatography techniques


- Hydrophobic interaction chromatography
 - Separates based on hydrophobicity
 - Elute using salt but this time
 - Start in high [salt]



Other chromatography techniques


- Hydrophobic interaction chromatography

- Separates based on hydrophobicity
- Elute using salt but this time
- Start in high [salt]
- Separate with decreasing [salt]



Other chromatography techniques




- Hydrophobic interaction chromatography
 - Choice of resins



Other chromatography techniques

- Hydrophobic interaction chromatography

- Choice of resins
- Increasing hydrophobicity for different separation



Chromatography Equipment

- Manual equipment – gravity columns
- Simple liquid handling - peristaltic pumps
- Automated systems – high performance liquid chromatography
 - Range of instruments from simple to complex
 - Range of suppliers
- All do essentially the same thing, pump liquid through a column



Manual equipment – gravity columns

- Best when using loose resin
- Manually load resin into columns
- Can run via gravity alone
 - You used these yesterday




Manual equipment – gravity columns

- Best when using loose resin
- Manually load resin into columns
- Can run via gravity alone
 - You used these yesterday
- Come in a range of sizes depending on the size of your prep and how much protein you have




Simple liquid handling – peristaltic pumps

- The theory




Simple liquid handling – peristaltic pumps

- The theory




Simple liquid handling – peristaltic pumps

- The theory



Simple liquid handling – peristaltic pumps

- The theory




Simple liquid handling – peristaltic pumps

- Come in a range of styles/models




Peristaltic pump for chromatography

- Can connect peristaltic pumps to:
 - gravity columns
 - small pre-packed columns




Peristaltic pump for chromatography

- Can connect peristaltic pumps to:
 - gravity columns
 - small pre-packed columns
- Need to:
 - manually collect elution
 - ensure sample/buffer doesn't run out




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components





Automated systems – AKTA Start

- Components





Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components




Automated systems – AKTA Start

- Components





Automated systems – AKTA Start

- Components



Automated systems – AKTA Start



- Sample Injection



Automated systems – AKTA Start

- Sample Injection



- Your protein sample will be in a syringe in position 3
- When the valve is set to LOAD, your sample can be loaded into the loop



Automated systems – AKTA Start

- Sample Injection

- Your protein sample will be in a syringe in position 3
- When the valve is set to LOAD, your sample can be loaded into the loop
- When the valve is set to INJECT, your sample (that is now in the loop) will be injected onto the column




Automated systems – AKTA Start


- Components
 - The system we'll be using has
 - a “carousel” fraction collector
 - computer control




Automated systems – AKTA Start



Automated systems – AKTA Pure



AKTA Chromatography Systems



Start



Go



Pure



Avant

Other Chromatography Systems



Agilent Technologies



Bio-Rad



Shimadzu



Pharmacia
(old AKTA)

Maintaining a chromatography system

- These systems and the associated columns are very sensitive to:
 - Particulates/aggregates



Maintaining a chromatography system

- These systems and the associated columns are very sensitive to:
 - Particulates/aggregates
 - Filter all buffers and samples



Maintaining a chromatography system

- These systems and the associated columns are very sensitive to:
 - Particulates/aggregates
 - Filter all buffers and samples
 - Air bubbles



Maintaining a chromatography system

- These systems and the associated columns are very sensitive to:
 - Particulates/aggregates
 - Filter all buffers and samples
 - Air bubbles
 - Degas buffers, keep bottles topped up
 - Careful protocols when injecting samples



Maintaining a chromatography system

- These systems and the associated columns are very sensitive to:
 - Particulates/aggregates
 - Filter all buffers and samples
 - Air bubbles
 - Degas buffers, keep bottles topped up
 - Careful protocols when injecting samples
 - High pressure



Maintaining a chromatography system

- These systems and the associated columns are very sensitive to:
 - Particulates/aggregates
 - Filter all buffers and samples
 - Air bubbles
 - Degas buffers, keep bottles topped up
 - Careful protocols when injecting samples
 - High pressure
 - Compression of the column resin is bad
 - Include pressure limit settings in protocols



Cleaning columns and safe storage

- After a chromatography run equilibrate your column back into:
 - Sterile-filtered (milli-Q) water, then
 - Ethanol, store in 20% ethanol
- Regularly clean the column with 0.2M NaOH
 - 0.2M NaOH will clean the column but is very corrosive so don't leave the column in this for any length of time
 - Either inject small volumes or use the Clean-In-Place (CIP) protocols
 - For this see your handbook or the relevant column datasheet



Protein Purification – Part 2

- Today we will learnt about:
 - Other chromatography techniques
 - Size-exclusion
 - Ion exchange
 - HIC and others
 - Chromatography equipment
 - The components and how to inject sample
 - The dos and don'ts of good chromatography



Questions?

Then let's do some SEC

